



## ACID DIGESTION

### CI ACID DIGESTION VESSEL

**CODE:08082**

The CI Acid Digestion Vessel (Bomb) helps in preparing analytical samples for analysis by combining various chemicals.

Inner Made of Teflon and outer made of Stainless steel with Brass plated lid.

### APPLICATIONS :

It is used in various industries like pharma, Pollution control boards, testing Labs, Research institutes, Universities. etc.. for digestion of Organic or Inorganic samples.



- \* Accelerated Digestion, as compared to open system at atmospheric pressure.
- \* Dissolves samples without losing original trace elements
- \* No unwanted Contaminations.
- \* Inner made of high density Teflon.
- \* Permits the use of strong acids, which will attack most other conventional containers.

MODEL	AD-23	AD-50	AD-100	AD-150
CODE	08082-00	08082-01	08082-02	08082-03
SIZE (ml)	23	50	100	150
SAMPLE VOL.(ml)	15	30	60	100
OUTER MATERIAL	STAINLESS STEEL			
INNER MATERIAL	TEFLON.			
MAX. TEMPERATURE	200°C			
CLOSURE	SPANNER		ALLEN SCREW (6)	

### ACCESSORIES/SPARES:

- \* Vice and Spanner Set
- \* Teflon inner Beaker
- \* Rupture Disc